

CLEANLINZ

CONCENTRATE

LIQUID CONCENTRATE

BEFORE



AFTER

Dental Evacuation Line Cleaning System
Restores System Suction to Full Capacity

BEFORE



AFTER

- Biodegradable
- Environmentally safe
- Eliminates sludge and reduces odor
- Non-Foaming and safe for dry pump systems
- Reduces the population of H2S forming bacteria

- Safe for all pipes
- All natural
- Safe for all animals & plants
- Provides optimized de-gradation of organic waste

Super concentrated formula digests organic matter & eliminates odors while improving vacuum suction.

CLEANLINZ breaks down and digests organic waste through the entire evacuation system which restores vacuum suction to its maximum capacity. The proprietary formula of bacteria colonizes through the evacuation system. A single cell will divide and multiply into 17 million cells in just 8 hours and continue to do this for approximately 8 days, spreading throughout the system and continuing to clean between uses.

Years of trapped accumulation left behind by other cleaners is released when **CLEANLINZ** is first used on the evacuation line system.

The first week of using **CLEANLINZ**, there may be increased sludge and debris accumulation in traps. This is the direct result of **CLEANLINZ** releasing the accumulation and shows the process of cleaning out the lines inner surfaces which will begin to restore suction. Therefore, the traps will have to be cleaned by hand until the system has been cleared of sludge. After using **CLEANLINZ** (depending on the age of the system), the traps will no longer require cleaning by hand, just dump traps once or twice per month.

100% Satisfaction Guaranteed

MARS
BIO-MED PROCESSES INC



facebook.com/MARSBIOEMED

575 Wentworth Street East, Unit 21, Oshawa, Ontario L1H 3V8
Tel: 905-723-9600 • 1-866-594-3648 • www.marsbiomed.com