

Line Cleaner on the pH Scale

Concentration of Hydrogen Ions compared to distilled water

10,000,000	pH= 0	PureVac (0.5-1.5)
1,000,000	pH= 1	BM 6000(1.05)
100,000	pH= 2	PureGreen 24 (2.0) AquaFx (2.5)
10,000	pH= 3	
1,000	pH= 4	
100	pH= 5	
10	pH= 6	Denta-Zyme (6.0-8.5)
1	pH= 7	CleanLINZ (6.6-7.4) PowerLINZ(7.0)
1/10	pH= 8	Neutravac (7.0-8.0) Em Power (7.5-8.5) Bio Vac(8.5)
1/100	pH= 9	CleanStream (8.5-10.9) PureVac Sc (9.4)
1/1,000	pH=10	MicroVac (9.0-10.0) VacAttak(10.5-11.5)
1/10,000	pH=11	Evaculine (11.5)
1/100,000	pH=12	Vacusol Ultra (12.4)
1/1,000,000		

preferred {