

# POWERLINZ<sup>®</sup>

## CONCENTRATE

### Evacuation Line Cleanser



M.A.R.S. Bio-Med Processes Inc. has created it's own brand of evacuation line cleaner to assure the maximum efficiency out of your amalgam separator. All evacuation line cleaners that are used with amalgam separators should have a neutral pH. Highly alkaline or acidic cleaners will release mercury and decrease the efficiency of the separator. In time, non-neutral cleaners will also negatively impact your evacuation lines and your vacuum pump system.

**PowerLINZ** is fast and effective. Our formulation controls the rapid bacterial build-up inside the suction lines along with the amalgam separator and is especially created for your *LibertyBOSS*.

#### Everything you need in a line cleaner:

- o Concentrated
- o Cost-Effective
- o Easy-to-Use
- o Convenient Premeasured Pump
- o Non-Foaming
- o Non-Corrosive
- o Antimicrobial
- o 2 Bottles/Case
- o Neutral pH
- o Chlorine-free
- o 100% Biodegradable
- o 1 Bottle makes 125 Liters



575 Wentworth Street E., Unit 21  
905-723-9600

Oshawa, Ontario  
[www.marsbiomed.com](http://www.marsbiomed.com)

L1H 3V8  
1-866-594-3648